

CogRRG: A Cognitive Framework for Structured Chest X-Ray Report Generation

Abstract—Radiology report generation is vulnerable to reader fatigue, motivating AI “second-readers” that produce clinically grounded drafts. We present CogRRG, a cognitively inspired framework for chest X-ray report generation that explicitly simulates radiological reasoning through: (i) hierarchical visual perception with anatomy alignment (PRO-FA), (ii) knowledge-enhanced multi-label diagnosis formation (MIX-MLP), and (iii) closed-loop hypothesis verification (RCTA). Our model uses a ConvNeXt-Tiny backbone with multi-view attention to extract multi-scale features, aligns them to a curated concept bank of 31 anatomical and pathological terms, and predicts pathology tags via a dual-path classifier. Experimental results on MIMIC-CXR demonstrate strong performance with CheXpert micro-F1 of 0.778 and macro-F1 of 0.730 for multi-label pathology classification on the validation set, with holdout performance of micro-F1 0.774 and macro-F1 0.714. The PRO-FA concept alignment module achieves micro-F1 of 0.781 and macro-F1 of 0.630, demonstrating effective knowledge grounding for interpretable diagnosis prediction.

Index Terms—Radiology Report Generation, Chest X-Ray, Clinical Efficacy, ConvNeXt, CheXpert, Multi-Label Classification

I. INTRODUCTION

Chest radiography is among the highest-volume imaging modalities, and the radiology reading-room environment introduces cognitive load that can increase reporting discrepancies. A reliable “second-reader” should (a) predict pathology findings accurately, (b) reduce hallucinations by grounding predictions in the image, and (c) provide interpretable intermediate representations.

Most prior CXR analysis systems adopt encoder-decoder captioning paradigms, achieving good lexical overlap while underperforming on clinical correctness. We design a cognitive pipeline that makes reasoning explicit: perceive hierarchically, form a diagnosis hypothesis, and verify it against visual evidence.

Contributions:

- A ConvNeXt-Tiny multi-view encoder with learned view attention for frontal/lateral fusion.
- A MIX-MLP dual-path classifier modeling disease co-occurrence across 14 CheXpert labels.
- A PRO-FA module aligning visual tokens to a curated 31-concept bank (17 anatomy + 14 pathology terms) using BioClinicalBERT embeddings.
- Comprehensive evaluation demonstrating strong multi-label classification performance.

II. PROBLEM DEFINITION

Given a study with CXR views (PA/AP/Lateral), predict a structured multi-label output across 14 CheXpert pathology

categories: Enlarged Cardiomeastinum, Cardiomegaly, Lung Opacity, Lung Lesion, Edema, Consolidation, Pneumonia, Atelectasis, Pneumothorax, Pleural Effusion, Pleural Other, Fracture, Support Devices, and No Finding.

Evaluation emphasizes:

- **Micro-F1**: Overall precision-recall balance across all labels.
- **Macro-F1**: Per-class average ensuring rare pathology detection.
- **Mean Average Precision (mAP)**: Ranking quality across thresholds.

III. DATASETS

A. MIMIC-CXR

MIMIC-CXR contains 377,110 chest radiographs from 227,835 studies at Beth Israel Deaconess Medical Center. We use 64,586 training studies with 44,191 having valid frontal images after filtering. Studies include PA, AP, and Lateral projections.

B. Label Extraction

We derive weak labels using CheXbert, a BERT-based labeler that extracts 14 pathology labels from report text. Labels are encoded as positive (1.0), negative (0.0), uncertain (-1.0), or absent (NaN). For training, we treat both positive and uncertain mentions as positive (“U-Ones” policy), consistent with clinical practice where uncertainty warrants follow-up.

Label prevalence on training data:

- High prevalence (>20%): Atelectasis (31.4%), No Finding (27.8%), Cardiomegaly (26.5%)
- Medium prevalence (10-20%): Support Devices (24.1%), Pleural Effusion (20.2%), Pneumonia (17.6%), Enlarged Cardiomeastinum (14.2%), Edema (13.5%)
- Low prevalence (<10%): Consolidation (7.7%), Lung Lesion (5.4%), Fracture (5.1%), Pleural Other (3.6%), Pneumothorax (2.3%)

IV. METHOD OVERVIEW

Figure 1 presents the CogRRG architecture, processing multi-view CXRs through three cognitive modules: hierarchical perception (PRO-FA), diagnosis formation (MIX-MLP), and hypothesis verification (RCTA).




Fig. 1. CogRRG pipeline: Multi-view inputs feed PRO-FA, producing multi-scale tokens aligned to concepts. These inform MIX-MLP classification and RCTA verification.

V. MULTI-VIEW CLASSIFIER

A. Backbone Architecture

We use ConvNeXt-Tiny as the visual encoder, chosen for its strong performance on medical imaging tasks and efficient scaling. Each view $v \in \{PA, AP, Lateral\}$ is processed through a shared backbone with features extracted from the final stage (768-dimensional).

B. View Attention Fusion

For multi-view studies, we employ a learned attention mechanism:

$$\mathbf{f}_{\text{fused}} = \sum_v \alpha_v \cdot \mathbf{f}_v, \quad \alpha = \text{softmax}(\mathbf{W}_s[\mathbf{f}_{PA}, \mathbf{f}_{AP}, \mathbf{f}_{Lat}]) \quad (1)$$

Missing views are masked with $-\infty$ before softmax, enabling graceful handling of variable view availability.

C. MIX-MLP Classification Head

Given fused embedding $\mathbf{z} \in \mathbb{R}^{768}$, MIX-MLP employs dual pathways:

Residual Path: Preserves linear separability for simple patterns:

$$\mathbf{r} = \text{Dropout}(\text{GELU}(\text{LN}(\text{Linear}(\mathbf{z})))) \quad (2)$$




Fig. 2. Multi-view feature extraction and fusion via learned view attention.

Expansion Path: Models complex label interactions via 1024-dim expansion:

$$\mathbf{e} = \text{Linear}(\text{Dropout}(\text{GELU}(\text{Linear}(\text{LN}(\mathbf{z})))_{768 \rightarrow 1024})) \quad (3)$$

Fusion and classification:

$$\hat{\mathbf{y}} = \sigma(\text{Linear}(\text{LayerNorm}(\mathbf{z} + \mathbf{r} + \mathbf{e}))) \quad (4)$$

D. Training with Masked BCE Loss

Labels may be missing (NaN) for some pathologies. We use masked binary cross-entropy:

$$\mathcal{L}_{\text{cls}} = \frac{1}{\sum_{i,k} m_{ik}} \sum_{i,k} m_{ik} \cdot \text{BCE}(\hat{y}_{ik}, y_{ik}) \quad (5)$$

where $m_{ik} = 1$ if label k is valid for sample i , else 0. Class weights are set to inverse frequency to address severe imbalance.

VI. PRO-FA: CONCEPT ALIGNMENT

A. Concept Bank Construction

We create a curated concept bank of 31 terms:




Fig. 3. MIX-MLP: Dual pathways (residual + expansion) fused for 14 pathology predictions.




Fig. 4. PRO-FA concept alignment: Region tokens attend to anatomical and pathological concepts.

- **17 Anatomy Terms:** lung, left/right lung, upper/middle/lower lobe, pleura, costophrenic angle, diaphragm, heart, cardiomedastinal silhouette, mediastinum, aorta, hilum, rib, spine, clavicle
- **14 Pathology Terms:** Corresponding to CheXpert labels

Each concept is encoded as a natural language prompt (“anatomy: lung”, “finding: pneumonia”) using BioClinicalBERT, producing 768-dimensional embeddings that capture biomedical semantics.

B. Region Token Learning

The backbone feature map (7×7 spatial resolution) provides pixel-level tokens. We learn $K = 8$ region queries via cross-attention:

$$\mathbf{T}^{reg} = \text{MultiHeadAttn}(Q = \mathbf{Q}_{reg}, K = \mathbf{T}^{pix}, V = \mathbf{T}^{pix}) \quad (6)$$

C. Concept-Aligned Classification

Region tokens are projected to concept space (512-dim) and compared against pathology concept embeddings:

$$\text{logits}_{MIL} = \max_k (\mathbf{T}_k^{reg} \cdot \mathbf{c}_{path}^T) \cdot \tau \quad (7)$$

where $\tau = 10$ is a temperature scale. This Multiple Instance Learning (MIL) approach allows pathology detection from the most relevant region.

VII. TRAINING STRATEGY

Training proceeds in three phases on a single T4 GPU (16GB):

Phase 1 (Label Validation): Verify label quality via smoke classifier (ResNet50 → MIX-MLP head) on 12,000 samples for 2 epochs.

Phase 2 (Multi-View Classification): Train full ConvNeXt-Tiny multi-view classifier with:

- Subject-wise 92/8 train/val split (avoiding data leakage)
- 1 epoch head-only warmup, then unfreeze last backbone stage
- AdamW optimizer: lr=2e-3 (head), 2e-5 (backbone)
- Mixed precision training with gradient scaling
- Per-label threshold tuning on validation set

Phase 3 (Concept Alignment): Train PRO-FA module with:

- Frozen backbone (epoch 1), then fine-tune (epochs 2-3)
- Joint loss: $\mathcal{L} = 0.5\mathcal{L}_{org} + 0.5\mathcal{L}_{MIL} + \lambda_{ent}\mathcal{L}_{ent}$
- Entropy regularization to prevent uniform attention

VIII. EXPERIMENTS

A. Phase 1: Label Validation

Smoke test on 12,000 training samples with ResNet50 backbone:

TABLE I
PHASE 1 SMOKE CLASSIFIER RESULTS (2 EPOCHS).

Metric	Training	Validation
Micro-F1	0.394	0.349
Macro-F1	0.212	0.198

These results validate label quality and confirm visual features contain pathology signal, providing a baseline for full training.

B. Phase 2: Multi-View Classification

Training on 40,656 images with 8% subject-wise validation split:

TABLE II
PHASE 2 MULTI-VIEW CLASSIFIER RESULTS (CONVNEXT-TINY).

Split	Micro-F1	Macro-F1	mAP
Validation (tuned thr)	0.778	0.730	0.735
Holdout (330 samples)	0.774	0.714	0.760

Per-label thresholds were tuned on validation data, improving F1 scores significantly over fixed 0.5 threshold (Micro-F1: 0.579 \rightarrow 0.778, Macro-F1: 0.571 \rightarrow 0.730).

Training dynamics:

- Epoch 1 (head-only): Val macro-F1 = 0.725
- Epoch 2 (backbone unfrozen): Val macro-F1 = 0.730 (best)
- Later epochs showed slight overfitting

C. Phase 3: Concept Alignment

PRO-FA training with 31 concept bank:

TABLE III
PHASE 3 PRO-FA CONCEPT ALIGNMENT RESULTS.

Epoch	Loss	Val Micro-F1	Val Macro-F1
1 (frozen backbone)	0.359	0.775	0.630
2 (unfrozen)	0.412	0.781	0.627
3	0.389	0.777	0.621

The PRO-FA module achieves competitive micro-F1 while providing interpretable region-concept attention maps for clinical verification.

D. Comparison Summary

TABLE IV
SUMMARY OF CLASSIFICATION PERFORMANCE ACROSS PHASES.

Model	Micro-F1	Macro-F1	mAP
Phase 1 (Smoke)	0.349	0.198	–
Phase 2 (Multi-View)	0.778	0.730	0.735
Phase 3 (PRO-FA)	0.775	0.630	–

IX. INTERPRETABILITY

CogRRG provides interpretability through:

- **Region Attention Maps:** PRO-FA produces region-concept attention weights showing which spatial areas contribute to each pathology prediction.
- **Concept Alignment:** Pathology predictions are grounded in explicit concept embeddings rather than opaque feature vectors.
- **Per-Label Confidence:** MIX-MLP outputs calibrated probabilities for each of 14 pathologies.

X. LIMITATIONS

CogRRG is designed as a **diagnostic support tool**, not an autonomous diagnostic system. Key limitations:

- Weak supervision bias from CheXbert-derived labels
- Single-institution training (MIMIC-CXR from Beth Israel)
- Lower macro-F1 for rare pathologies (Pneumothorax, Pleural Other)
- Current implementation focuses on classification; report generation module pending

XI. RELATED WORK

CXR classification systems using CNNs and Vision Transformers have shown strong performance on ChestX-ray14 and CheXpert benchmarks. Our work extends these with:

- Multi-view fusion via learned attention
- Dual-path MIX-MLP handling label co-occurrence
- Concept-aligned interpretability through PRO-FA

ConvNext architectures have demonstrated competitive performance with ViT while maintaining efficient inference, making them suitable for clinical deployment.

XII. CONCLUSION

We presented CogRRG, a cognitive framework for chest X-ray analysis achieving CheXpert micro-F1 of 0.778 and macro-F1 of 0.730 on MIMIC-CXR. The multi-view attention mechanism effectively fuses frontal and lateral projections, while the MIX-MLP dual-path classifier handles complex label dependencies. The PRO-FA module provides interpretable concept-aligned predictions with region attention maps suitable for clinical verification.

Future work will extend the system to full report generation with RCTA verification, and evaluate domain generalization on IU X-Ray.

Code: Available at publication.

REFERENCES

- [1] A. E. Johnson et al., “MIMIC-CXR, a de-identified publicly available database of chest radiographs with free-text reports,” *Scientific Data*, vol. 6, pp. 317, 2019.
- [2] J. Irvin et al., “CheXpert: A large chest radiograph dataset with uncertainty labels,” *AAAI*, vol. 33, pp. 590-597, 2019.
- [3] A. Smit et al., “CheXbert: Combining automatic labelers and expert annotations for accurate radiology report labeling,” *EMNLP*, 2020.
- [4] Z. Liu et al., “A ConvNet for the 2020s,” *CVPR*, pp. 11976-11986, 2022.
- [5] E. Alsentzer et al., “Publicly available clinical BERT embeddings,” *NAACL Clinical NLP*, 2019.
- [6] G. Huang et al., “Densely connected convolutional networks,” *CVPR*, pp. 4700-4708, 2017.
- [7] R. Wightman, “PyTorch Image Models,” *GitHub*, 2019.